

2021/09



NEW PRODUCTS



LIPOMED continuously produces new reference standards to assist our customers with the ever increasing need for high quality and reliable reference standards.


We exclusively handle the full production and release process of new products to ensure highest quality. Synthesis, stability testing, quality control as well as compounding, labeling and distribution are all performed at our headquarters in Switzerland.

Should an item you are looking for not be included in our current portfolio please contact our customer service (see last page) or visit our website for more information.

- » Alcohol Mixture 4 (MIX-1808-2LH)
- » Alcohol Mixture 4 (MIX-1808-5LH)
- » Clozapine-D₄
- » Dehydroaripiprazole
- » N-Desmethyloanzapine
- » 5F-APINACA
- » 3-(4-Fluorobenzoyl)-propionic acid
- » Hydroxybupropion-D₆
- » Iso-LSD
- » Mesoridazine.besylate
- » Nortriptyline-D₃.HCl
- » 2-Oxo-3-hydroxy-LSD
- » Protriptyline-D₃.HCl
- » Trazodone-D₆.HCl
- » Vilazodone-D₄.HCl
- » Vortioxetine-D₈.HBr


new

<p>Acetone Ethanol Isopropanol Methanol</p>	<p>Alcohol Mixture 4</p>		
	<p>MIX-1808-2LH</p>	<p>2.0 mg/1 ml of each component</p>	<p>calibrated solution in water</p>
<p>Acetone Ethanol Isopropanol Methanol</p>	<p>Alcohol Mixture 4</p>		
	<p>MIX-1808-5LH</p>	<p>5.0 mg/1 ml of each component</p>	<p>calibrated solution in water</p>
	<p>Clozapine-D₄</p>		
	<p>C₁₈H₁₅D₄N₄</p>	<p>M.W.: 330.85</p>	
	<p>CLZ-1913-FB-0.1LM</p>	<p>0.1 mg base/1 ml</p>	<p>calibrated solution in methanol</p>
	<p>CLZ-1913-FB-1LM</p>	<p>1.0 mg base/1 ml</p>	<p>calibrated solution in methanol</p>
	<p>Dehydroaripiprazole</p>		
	<p>C₂₃H₂₅Cl₂N₃O₂</p>	<p>M.W.: 446.37</p>	
	<p>ARP-1883-1LM</p>	<p>1.0 mg/1 ml</p>	<p>calibrated solution in methanol</p>
	<p>ARP-1883-10</p>	<p>10 mg</p>	<p>powder</p>
	<p>ARP-1883-50</p>	<p>50 mg</p>	<p>powder</p>
	<p>ARP-1883-100</p>	<p>100 mg</p>	<p>powder</p>




N-Desmethylolanzapine

$C_{16}H_{18}N_4S$	M.W.: 298.41	
OLZ-1899-FB-1LA	1.0 mg base/1 ml	calibrated solution in acetonitrile
OLZ-1899-FB-10	10 mg	powder
OLZ-1899-FB-50	50 mg	powder
OLZ-1899-FB-100	100 mg	powder




5F-APINACA

Synonym(s): 5F-AKB-48		
$C_{23}H_{30}FN_3O$	M.W.: 383.51	
AKB-1910-0.1LM	0.1 mg/1 ml	calibrated solution in methanol
AKB-1910-1LM	1.0 mg/1 ml	calibrated solution in methanol
AKB-1910-10	10 mg	powder
AKB-1910-50	50 mg	powder
AKB-1910-100	100 mg	powder




3-(4-Fluorobenzoyl)-propionic acid

Synonym(s): Haloperidol metabolite III		
$C_{10}H_9FO_3$	M.W.: 196.18	
FPA-1879-FA-1LM	1.0 mg acid/1 ml	calibrated solution in methanol
FPA-1879-FA-10	10 mg	powder
FPA-1879-FA-50	50 mg	powder
FPA-1879-FA-100	100 mg	powder



Hydroxybupropion-D₆

$C_{13}H_{12}D_6ClNO_2$	M.W.: 261.78	
HBP-1912-FB-0.1LA	0.1 mg base/1 ml	calibrated solution in acetonitrile
HBP-1912-FB-1LA	1.0 mg base/1 ml	calibrated solution in acetonitrile
HBP-1912-FB-10	10 mg	powder
HBP-1912-FB-50	50 mg	powder
HBP-1912-FB-100	100 mg	powder




Iso-LSD

Synonym(s): Iso-Lysergic acid diethylamide

$C_{20}H_{25}N_3O$ M.W.: 323.44


LSD-1944-FB-0.1LA	0.1 mg base/1 ml	calibrated solution in acetonitrile
LSD-1944-FB-1LA	1.0 mg base/1 ml	calibrated solution in acetonitrile
LSD-1944-FB-10	10 mg	powder
LSD-1944-FB-50	50 mg	powder
LSD-1944-FB-100	100 mg	powder



Mesoridazine.besylate

$C_{21}H_{26}N_2O_2 \cdot C_6H_6O_3S$ M.W.: 386.57 + 158.17


MES-1880-BS-10	10 mg	powder
MES-1880-BS-50	50 mg	powder
MES-1880-BS-100	100 mg	powder



Nortriptyline-D₃.HCl

$C_{19}H_{18}D_3N.HCl$ M.W.: 266.40 + 36.46

NTP-1909-HC-0.1LM	0.1 mg base/1 ml	calibrated solution in methanol
NTP-1909-HC-1LM	1.0 mg base/1 ml	calibrated solution in methanol
NTP-1909-HC-10	10 mg	powder
NTP-1909-HC-50	50 mg	powder
NTP-1909-HC-100	100 mg	powder




2-Oxo-3-hydroxy-LSD

Synonym(s): 2-Oxo-3-hydroxy-lysergic acid diethylamide

$C_{20}H_{25}N_3O_3$ M.W.: 355.43


LSD-1946-0.1LA	0.1 mg/1 ml	calibrated solution in acetonitrile
LSD-1946-1LA	1.0 mg/1 ml	calibrated solution in acetonitrile
LSD-1946-10	10 mg	powder
LSD-1946-50	50 mg	powder
LSD-1946-100	100 mg	powder



Protriptyline-D₃.HCl

C₁₉H₁₈D₃N.HCl M.W.: 266.40 + 36.46


PTR-1901-HC-0.1LM	0.1 mg base/1 ml	calibrated solution in methanol
PTR-1901-HC-1LM	1.0 mg base/1 ml	calibrated solution in methanol
PTR-1901-HC-10	10 mg	powder
PTR-1901-HC-50	50 mg	powder
PTR-1901-HC-100	100 mg	powder



Trazodone-D₆.HCl

C₁₉H₁₆D₆ClN₅O.HCl M.W.: 377.90 + 36.46


TZD-1904-HC-0.1LM	0.1 mg base/1 ml	calibrated solution in methanol
TZD-1904-HC-1LM	1.0 mg base/1 ml	calibrated solution in methanol
TZD-1904-HC-10	10 mg	powder
TZD-1904-HC-50	50 mg	powder
TZD-1904-HC-100	100 mg	powder



Vilazodone-D₄.HCl

C₂₆H₂₃D₄N₅O₂.HCl M.W.: 445.56 + 36.46

VZD-1905-HC-0.1LAH	0.1 mg base/1 ml	cal. solution in acetonitrile/water (1:1)
VZD-1905-HC-1LAH	1.0 mg base/1 ml	cal. solution in acetonitrile/water (1:1)
VZD-1905-HC-10	10 mg	powder
VZD-1905-HC-50	50 mg	powder
VZD-1905-HC-100	100 mg	powder



Vortioxetine-D₈.HBr

C₁₈H₁₈D₈N₂S.HBr M.W.: 306.50 + 80.91

VOX-1916-HB-0.1LM	0.1 mg base/1 ml	calibrated solution in methanol
VOX-1916-HB-1LM	1.0 mg base/1 ml	calibrated solution in methanol
VOX-1916-HB-10	10 mg	powder
VOX-1916-HB-50	50 mg	powder
VOX-1916-HB-100	100 mg	powder

LIPOMED MAIL/FAX Quotation & Order Form

Contact Information

Bill to	Ship to (if different from invoicing address)
Customer/Company	Customer/Company
Department	Department
Address	Address
Contact Person	Contact Person
City, ZIP	City, ZIP
State, Country	State, Country
Phone	Phone
Fax	Fax
E-Mail Address	E-Mail Address
No of pages (incl. this one)	Customer Account #
Date of Request	Purchase Order
Special instructions/Notes	

Product List

Article-No.	Product Name	Qty	Unit Price	Total Price

Fax France +33 3 88 96 50 73 | Germany +49 7621 1693 474 | USA +1 617 577 1776 | CH & other countries +41 61 702 02 20
E-Mail USA ussales@lipomed.com | All other countries lipomed@lipomed.com
In case a permit is needed, please attach original permit along with your order and send it to Lipomed per mail.
Reference Standards orders are subject to LIPOMED's standard Terms and Conditions of Sale.

Ordering & other information

LIPOMED is a private Swiss Health Care Company. We are one of the leading companies for Analytical Reference Standards worldwide and specialized in the research, development and manufacturing of ethical pharmaceuticals and products for the treatment of rare diseases.

Since 1993 LIPOMED is dedicated to providing products and services of outstanding quality and expertise. We have our own GMP/GDP and ISO 9001 certified pharmaceutical production facilities in Arlesheim, near Basel. We are also ISO/IEC 17025 accredited for testing of Analytical Reference Standards and ISO 17034 accredited for the production of Reference Material.

Our innovative products are marketed and distributed with offices and partners in more than 50 countries around the world.

LIPOMED's Analytical Reference Standards currently include the following classes of substances, their deuterated analogues and metabolites, such as Glucuronides and Sulfates:

- » Alcohol & Ethanol Biomarkers
- » Amphetamines & Sympathomimetic amines
- » Antidepressants
- » Antipsychotics
- » Barbiturates
- » Benzodiazepines
- » Cannabinoids
- » Cathinones
- » Cocaines
- » Custom Mixtures
- » Hallucinogens
- » Opioids
- » Other CNS Depressants
- » Steroids

For more information about LIPOMED's Analytical Reference Standards, please visit our website www.lipomed.com.

USA, CANADA, BERMUDA

LIPOMED Inc
150 Cambridgepark Drive
Cambridge, MA 02140
USA

Tel +1 617 577 7222
Fax +1 617 577 1776
ussales@lipomed.com

SWITZERLAND (HEADQUARTERS)

LIPOMED AG
Fabrikmattenweg 4
4144 Arlesheim
Switzerland

Tel +41 61 702 02 00
Fax +41 61 702 02 20
lipomed@lipomed.com

GERMANY

LIPOMED GmbH
Hegenheimer Str. 2
79576 Weil am Rhein
Germany

Tel +49 7621 1693 473
Fax +49 7621 1693 474
lipomed@lipomed.com



LIPOMED AG, Switzerland is **GMP/GDP** and **ISO 9001** certified for production and distribution of Pharmaceuticals and **ISO/IEC 17025** (ANAB certificate no. AT-1760) accredited for testing of Analytical Reference Standards and **ISO 17034** (ANAB certificate no. AR-1761) accredited for Reference Material production. All products distributed by LIPOMED affiliates are manufactured by LIPOMED AG and distributed under protocols of LIPOMED AG.

