



lipomed
SERVICES TO HEALTH

Analytical Reference Standards

2023/03



NEW PRODUCTS

new


LIPOMED continuously produces new reference standards to assist our customers with the ever increasing need for high quality and reliable reference standards.


We exclusively handle the full production and release process of new products to ensure highest quality. Synthesis, stability testing, quality control as well as compounding, labeling and distribution are all performed at our headquarters in Switzerland.


Should an item you are looking for not be included in our current portfolio please contact our customer service (see last page) or visit our website for more information.


- » Barbitol-D₅
- » Biperiden.HCl
- » Cannabigerolic acid
- » 1-(3-Chlorophenyl)-piperazine-D₈.HCl
- » Clozapine-D₄
- » Cyclobenzaprine-D₃.HCl
- » Digoxin
- » Ecgonine-D₃.HCl
- » Flualprazolam
- » Hydroxybupropion-D₆
- » alpha-Hydroxyflualprazolam
- » Iso-LSD
- » 5-MeO-DMT.HCl
- » (S)-Norfluoxetine.oxalate
- » 2-Oxo-3-hydroxy-LSD
- » d,l-PMA-D₃.HCl
- » Thiopental-D₅
- » Pentylone-D₃.HC
- » Valproic acid-D₆
- » Vortioxetine-D₈.HBr
- » Warfarin
- » Warfarin, 1.0 mg/ml solution
- » Zuclopenthixol


www.lipomed-usa.com


	Barbitol-D₅				
	Molecular Formula	C ₈ H ₇ D ₅ N ₂ O ₃	CAS Number	1189694-78-7	Status
	Molecular Weight	189.23	DEA Number	2145 CIV	CH DE FR
	BRB-1943-FA-0.1LM	0.1 mg acid/1 mL	calibrated solution in methanol		A D
	BRB-1943-FA-1LM	1.0 mg acid/1 mL	calibrated solution in methanol		A D
	BRB-1943-FA-10	10 mg	powder		A A D
	BRB-1943-FA-50	50 mg	powder		A A D
	BRB-1943-FA-100	100 mg	powder		A A D


	Biperiden.HCl				
	Molecular Formula	C ₂₁ H ₂₉ NO.HCl	CAS Number	1235-82-1	Status
	Molecular Weight	311.47 + 36.46	DEA Number		CH DE FR
	BIP-1918-HC-1LM	1.0 mg base/1 mL	calibrated solution in methanol		
	BIP-1918-HC-10	10 mg	powder		
	BIP-1918-HC-50	50 mg	powder		
	BIP-1918-HC-100	100 mg	powder		


	Cannabigerolic acid				
	Synonym(s):	CBGA	CAS Number	25555-57-1	Status
	Molecular Formula	C ₂₂ H ₃₂ O ₄	DEA Number	N/A	CH DE FR
	Molecular Weight	360.49			
	CBG-1924-10	10 mg	powder		
	CBG-1924-50	50 mg	powder		
	CBG-1924-100	100 mg	powder		


	1-(3-Chlorophenyl)-piperazine-D₈.HCl (shipment under ice for solution)				
	Synonym(s):	mCPP-D ₈ .HCl	CAS Number	1313393-63-3	Status
	Molecular Formula	C ₁₀ H ₅ D ₈ ClN ₂ .HCl	DEA Number	N/A	CH DE FR
	Molecular Weight	204.73 + 36.46			
	PPZ-1920-HC-0.1LM	0.1 mg base/1 mL	calibrated solution in methanol		A
	PPZ-1920-HC-1LM	1.0 mg base/1 mL	calibrated solution in methanol		A
	PPZ-1920-HC-10	10 mg	powder		A A
	PPZ-1920-HC-50	50 mg	powder		A A
	PPZ-1920-HC-100	100 mg	powder		A A


	Clozapine-D₄				
	Molecular Formula	C ₁₈ H ₁₃ D ₄ ClN ₄	CAS Number	N/A	Status
	Molecular Weight	330.85	DEA Number	N/A	CH DE FR
	CLZ-1913-FB-0.1LM	0.1 mg base/1 mL	calibrated solution in methanol		
	CLZ-1913-FB-1LM	1.0 mg base/1 mL	calibrated solution in methanol		


	Cyclobenzaprine-D₃.HCl			
	Molecular Formula	C ₂₀ H ₁₈ D ₃ N.HCl	CAS Number	1184983-42-3
	Molecular Weight	278.41 + 36.46	DEA Number	N/A
	CBZ-1926-HC-0.1LM	0.1 mg base/1 mL	calibrated solution in methanol	
	CBZ-1926-HC-1LM	1.0 mg base/1 mL	calibrated solution in methanol	
	CBZ-1926-HC-10	10 mg	powder	
	CBZ-1926-HC-50	50 mg	powder	
	CBZ-1926-HC-100	100 mg	powder	

	Digoxin			
	Molecular Formula	C ₄₁ H ₆₄ O ₁₄	CAS Number	20830-75-5
	Molecular Weight	780.94	DEA Number	N/A
	DOX-1921-1LM	1.0 mg base/1 mL	calibrated solution in methanol	
	DOX-1921-10	10 mg	powder	
	DOX-1921-50	50 mg	powder	
	DOX-1921-100	100 mg	powder	

	Ecgonine-D₃.HCl (shipment under ice for solution)			
	Molecular Formula	C ₉ H ₁₂ D ₃ NO ₃ .HCl	CAS Number	N/A
	Molecular Weight	188.24 + 36.46	DEA Number	9180 CII
	COC-295-HC-0.1LM	0.1 mg base/1 mL	calibrated solution in methanol	A C
	COC-295-HC-1LM	1.0 mg base/1 mL	calibrated solution in methanol	A C
	COC-295-HC-10	10 mg	powder	A A C
	COC-295-HC-50	50 mg	powder	A A C
	COC-295-HC-100	100 mg	powder	A A C


	Flualprazolam			
	Molecular Formula	C ₁₇ H ₁₂ ClFN ₄	CAS Number	28910-91-0
	Molecular Weight	326.76	DEA Number	N/A
	FAL-1919-FB-1LM	1.0 mg base/1 mL	calibrated solution in methanol	A A(*) D
	FAL-1919-FB-10	10 mg	powder	A A(*) D
	FAL-1919-FB-50	50 mg	powder	A A(*) D
	FAL-1919-FB-100	100 mg	powder	A A(*) D

	Hydroxybupropion-D₆			
	Molecular Formula	C ₁₃ H ₁₂ D ₆ ClNO ₂	CAS Number	N/A
	Molecular Weight	261.78	DEA Number	N/A
	HBP-1912-FB-0.1LA	0.1 mg base/1 mL	calibrated solution in acetonitrile	
	HBP-1912-FB-1LA	1.0 mg base/1 mL	calibrated solution in acetonitrile	
	HBP-1912-FB-10	10 mg	powder	
	HBP-1912-FB-50	50 mg	powder	
	HBP-1912-FB-100	100 mg	powder	



alpha-Hydroxyflualprazolam


Molecular Formula	C ₁₇ H ₁₂ ClFN ₄ O	CAS Number	37115-46-1	Status
Molecular Weight	342.76	DEA Number	N/A	CH DE FR
FAL-1925-FB-1LM	1.0 mg/1 mL	calibrated solution in methanol		
FAL-1925-FB-10	10 mg	powder		
FAL-1925-FB-50	50 mg	powder		
FAL-1925-FB-100	100 mg	powder		



Iso-LSD

Synonym(s): Iso-Lysergic acid diethylamide


Molecular Formula	C ₂₀ H ₂₅ N ₃ O	CAS Number	N/A	Status
Molecular Weight	323.44	DEA Number	N/A	CH DE FR
LSD-1944-FB-0.1LA	0.1 mg base/1 mL	calibrated solution in acetonitrile		
LSD-1944-FB-1LA	1.0 mg base/1 mL	calibrated solution in acetonitrile		
LSD-1944-FB-10	10 mg	powder		
LSD-1944-FB-50	50 mg	powder		
LSD-1944-FB-100	100 mg	powder		



5-MeO-DMT.HCl


Synonym(s): 5-Methoxy-N,N-dimethyltryptamine.HCl

Molecular Formula	C ₁₃ H ₁₈ N ₂ O.HCl	CAS Number	2427-79-4	Status
Molecular Weight	254.76	DEA Number	7431 CI	CH DE FR
DMT-1957-HC-1LM	1.0 mg base/1 mL	calibrated solution in methanol		
DMT-1957-HC-10	10 mg	powder		
DMT-1957-HC-50	50 mg	powder		
DMT-1957-HC-100	100 mg	powder		



(S)-Norfluoxetine.oxalate

Molecular Formula	C ₁₆ H ₁₀ D ₆ F ₃ NO.C ₂ H ₂ O ₄	CAS Number	126924-38-7 (FB)	Status
Molecular Weight	295.31 + 90.03	DEA Number	N/A	CH DE FR
FLX-1963-OX-1LM	1.0 mg base/1 mL	calibrated solution in methanol		
FLX-1963-OX-10	10 mg	powder		
FLX-1963-OX-50	50 mg	powder		
FLX-1963-OX-100	100 mg	powder		




2-Oxo-3-hydroxy-LSD

Synonym(s): 2-Oxo-3-hydroxy-lysergic acid diethylamide

Molecular Formula	C ₂₀ H ₂₅ N ₃ O ₃	CAS Number	N/A	Status
Molecular Weight	355.43	DEA Number	N/A	CH DE FR
LSD-1946-0.1LA	0.1 mg/1 mL	calibrated solution in acetonitrile		
LSD-1946-1LA	1.0 mg/1 mL	calibrated solution in acetonitrile		
LSD-1946-10	10 mg	powder		
LSD-1946-50	50 mg	powder		
LSD-1946-100	100 mg	powder		

LGC worldwide



1 Brazil
T +55 12 3302 5880
E bz@lgcstandards.com

2 Bulgaria
T +359 (0)2 971 4955
E bg@lgcgroup.com

3 China
T +86 400 9216156
E info.china@lgcgroup.com

4 France
T +33 (0)3 88 04 82 82
E fr@lgcstandards.com

5 Germany
T +49 (0)281 9887 0
E de@lgcstandards.com

6 Hungary
T +49 (0)281 9887 0
E de@lgcstandards.com

7 India
T +91 9082974025
E india@lgcgroup.com

8 Ireland
T +44 (0)20 8943 8480
E uksales@lgcstandards.com

9 Italy
T +39 02 22476412
E it@lgcstandards.com

10 Middle East
T +49 (0) 281 9887 270
E global.sales@lgcgroup.com

11 Netherlands
T +49 (0)281 9887 250
E de@lgcstandards.com

12 Nordic Countries
T +49 (0)281 9887 0
E de@lgcstandards.com

13 Poland
T +48 22 751 31 40
E pl@lgcgroup.com

14 Romania
T +49 (0)281 9887 0
E de@lgcstandards.com

15 Russia
T +7 812 777 04 88
E ru@lgcgroup.com

16 South Africa
T +27 (0)11 466 4321
E sales.za@lgcgroup.com

17 Spain
T +34 (0)93 308 4181
E es@lgcstandards.com

18 UK reference materials
T +44 (0)20 8943 8480
E uksales@lgcstandards.com

19 USA + Canada
T +1 (603) 622 7660
E lgcusa@lgcgroup.com

Export queries
T +49 (0) 281 9887 270
E global.sales@lgcgroup.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2023. All rights reserved.

Copyright © LGC Limited 2023. All rights reserved. LGC is a trademark owned by LGC GmbH